

ANTHEM[®]

Proximal Tibia Fracture System



Comprehensive Anatomic Plating Portfolio



XR Lateral Locking



Lateral Locking



Lateral Non-Locking



Medial Locking




Posteromedial Buttress

YOUR COMPLETE SOLUTION

Polyaxial design
±20° cone of angulation



Triple kickstand
Allows targeted and extended
medial column support



Dual rafting row
Provides additional rafting
support of articular surface

Radiolucent Instruments
Aid in fracture site visibility



GLOBUS
M E D I C A L